



## Digi-Flex®

Rehabilitative exercise and stress relief

### Characteristics

Function	The Digi-Flex® strengthens hands, fingers, forearms and will work the muscles separately or in tandem.
Application	<p>The Digi-Flex® can be used for rehab as well as for sports training or stress relief.</p> <ul style="list-style-type: none"> <li>• Rehabilitate post-operative / post-traumatic hands and arms</li> <li>• Restore function in hands, arms and shoulder following stroke</li> <li>• Reduce pain and stiffness in arthritic fingers, hands and wrists</li> <li>• Relax tense, sore muscles and reduce stress</li> <li>• Strengthen hands and arms</li> </ul>

### Technical specifications

Productgroup	Hand Therapy Products
Composition	Hard plastic body with cushioned rubber and spring loaded buttons

### Storage conditions

Duration of storage	No limitation
Storage	Store in cool, dry place.

## Product overview






Picture	Art.Nr	Description	Resistance	Package	Weight Package	Dimensions L x W x H
	DF150	Digi-Flex Yellow	0,7 kg / 2,3 kg	1 piece / sac	0,09 kg	12,5 x 8 x 3 cm
	DF300	Digi-Flex Red	1,4 kg / 4,5 kg	1 piece / sac	0,09 kg	12,5 x 8 x 3 cm
	DF500	Digi-Flex Green	2,3 kg / 7,3 kg	1 piece / sac	0,09 kg	12,5 x 8 x 3 cm
	DF700	Digi-Flex Blue	3,5 kg / 10,4 kg	1 piece / sac	0,09 kg	12,5 x 8 x 3 cm
	DF900	Digi-Flex Black	4,1 kg / 14,1 kg	1 piece / sac	0,09 kg	12,5 x 8 x 3 cm
	DFS10	Digi-Flex Set + display	5 resistances	1 set / box	1,20 kg	21 x 19 x 11 cm
	DF001	Empty Display		1 piece / box	1,00 kg	21 x 19 x 11 cm


## Application

### Safety Precautions

Use the Digi-Flex® only after consulting a trained licensed healthcare professional.

## Exercises



Special items	
Customs code	9021.9090
Country of origin	Korea MSD Europe is the authorised European Representative.
	This product has been manufactured according to CE-mark regulations.
Revision number	PI-Digi-Flex (04.2008)
<b>MSD Europe bvba</b> Neringstraat 7 1840 Londerzeel Belgium Tel. +32 52 311 452 Fax. +32 52 311 450 info@msd-europe.com www.msd-europe.com	