

MSD Hot/Cold Packs

Hot / Cold Therapy Treatment



Characteristics

Function

The MSD Hot/Cold packs are usable both when treating with heat and/or coldness. Therefore it is a strong alternative to the other ordinary heating pillows and ice packs.

The versatile, re-usable gel packs deliver hot and cold therapy effectively and economically. Even after being frozen under the 20° C it remains soft.

The MSD Hot/Cold Packs can be bent for any angle for the joint and muscle.

Technical specifications


Productgroup


Hot/Cold Therapy Products

Composition

The special gel is non-toxic and biodegradable

Product overview

Picture	Art.Nr	Description	Package	Weight Package	Dimensions L x W x H
	MSD-200	Small (15 x 25 cm)	10 pieces / box	6,9 kg	25 x 27 x 13 cm
	MSD-400	Medium (35 x 25 cm)	10 pieces / box	10,0 kg	35 x 24 x 17 cm
	MSD-600	Tri-Sectional (40 x 20 cm)	10 pieces / box	10,6 kg	41 x 24 x 16 cm
	MSD-900	Large (50 x 20 cm)	10 pieces / box	12,8 kg	50 x 20 x 16 cm

Application	
Safety Precautions	Use the MSD Hot/Cold Pack only after consulting a trained licensed healthcare professional.
Storage conditions	
Duration of storage	Valid for 3 years from date of manufacture
Storage	Store in cool, dry place, avoid excessive heat (40° C.)
Cleaning	They are reusable, washable and are very easy to disinfect
Special items	
Customs code	3005.9099
Country of origin	China MSD Europe is the authorised European Representative.
	This product has been manufactured according to CE-mark regulations.
Revision number	PI-MSD Hot-Cold Packs (04.2008)
MSD Europe bvba Neringstraat 7 1840 Londerzeel Belgium Tel. +32 52 311 452 Fax. +32 52 311 450 info@msd-europe.com www.msd-europe.com	