



MSD Instant Hot Pack

Instant Hot Therapy Treatment

Characteristics

Function	<p>The MSD Instant Hot Pack is a disposable pack, one-time use only. This pack is used for immediate hot therapy.</p> <p>By finger pressure the pack will activate. The insulation on the body side of the pack prevents dissipation of heat. The therapeutic pack maintains its temperature for about 30 minutes and the approximate peak temperature is 43° C.</p>
----------	--

Technical specifications

Productgroup	Hot / Cold Therapy Products
Composition	Magnesium Sulphate and Water

Product overview


Art.Nr	Description	Package	Weight Package	Dimensions L x W x H
MSD-250	MSD Instant Hot Pack	10 pieces / box	2,7 kg	18 x 24 x 15 cm

Storage conditions

Duration of storage	Valid for 3 years from date of manufacture
Storage	Store in cool, dry place, avoid excessive heat (40° C.)
Cleaning	Rince easily with water and soap

Application

Safety Precautions	Use the MSD Instant Hot Pack only after consulting a trained licensed healthcare professional.
--------------------	--

Special items	
Customs code	3005.9099
Country of origin	China MSD Europe is the authorised European Representative.
	This product has been manufactured according to CE-mark regulations.
Revision number	PI-MSD Instant Hot Pack (04.2008)
MSD Europe bvba Neringstraat 7 1840 Londerzeel Belgium Tel. +32 52 311 452 Fax. +32 52 311 450 info@msd-europe.com www.msd-europe.com	