



MSD Soft Touch Hot/Cold Packs

Hot/Cold Therapy Products

NO WRAP - NO COVER


Characteristics

Function	<p>The MSD Soft Touch Hot/Cold packs are usable when treating with both heat and/or coldness. Therefore it is a strong alternative to the other ordinary heating pillows and ice packs.</p> <p>The versatile, re-usable gel packs deliver hot and cold therapy effectively and economically. Even after being frozen under the 20° C it remains soft. The Packs can be used directly to the skin, NO COVER, NO WRAP.</p> <p>The MSD Hot/Cold Packs can be bent for any angle for the joint and muscle.</p>
Application	<ul style="list-style-type: none"> • (Sport) injuries • Tendinitis • Swelling • Muscle soreness • Strains • Sprains

Technical specifications

Productgroup	Hot/Cold Therapy Products
Composition	The special gel is non-toxic and biodegradable

Product overview

Picture	Art.Nr	Description	Package	Weight Package	Dimensions L x W x H
	MSD-650	Small (15 x 25 cm)	10 pieces / box	6,9 kg	25 x 27 x 13 cm
	MSD-950	Medium (35 x 25 cm)	10 pieces / box	10,0 kg	35 x 24 x 17 cm
	MSD-1000	Tri-Sectional (40 x 20 cm)	10 pieces / box	10,7 kg	41 x 24 x 16 cm
	MSD-1200	Large (50 x 20 cm)	10 pieces / box	12,5 kg	50 x 20 x 16 cm

Application

Safety Precautions Use the MSD Soft Touch Hot/Cold Pack only after consulting a trained licensed healthcare professional.

Storage conditions

Duration of storage Valid for 3 years from date of manufacture


Storage Store in cool, dry place, avoid excessive heat (40° C.)

Cleaning They are reusable, washable and are very easy to disinfect

Special items

Customs code 3005.9099

Country of origin China
MSD Europe is the authorised European Representative.

 This product has been manufactured according to CE-mark regulations.

Revision number PI-MSD Soft Touch Hot-Cold Packs (04.2008)

MSD Europe bvba
Neringstraat 7
1840 Londerzeel
Belgium
Tel. +32 52 311 452
Fax. +32 52 311 450
info@msd-europe.com
www.msd-europe.com