

Sit & Reach Unit

Evaluation and Measuring Device




Characteristics

Function	When using the Sit & Reach unit you sit with your feet against the “sit and reach box” and flex forward moving the slide bar. Flexion is read directly from the calibrated bar.
----------	---

Technical specifications

Productgroup	Evaluation line Instruments
Composition	Metal

Product overview


Picture	Art.Nr	Description	Package	Weigt Package	Dimensions L x W x H
	12-1085	Sit & Reach Unit	1 piece / case	5,00 kg	52.5 x 31 x 4.5 cm

Storage conditions

Duration of storage	No limitation
Storage	No special requirements
Cleaning	Rince easily with water and soap

Application

Safety Precautions	Use the Sit & Reach Unit only after consulting a trained licensed healthcare professional.
--------------------	--

Special items	
Customs code	9031.8098
Country of origin	U.S.A. MSD Europe is the authorised European Representative.
	This product has been manufactured according to CE-mark regulations.
Revision number	PI-Sit & reach Unit (04.2008)
MSD Europe bvba Neringstraat 7 1840 Londerzeel Belgium Tel. +32 52 311 452 Fax. +32 52 311 450 info@msd-europe.com www.msd-europe.com	