



## Theraflex™ Putty

Rehabilitative exercise and stress relief

### Characteristics

Function	The Theraflex™ Putty strengthens weakened grasp through exercise, is used for rehabilitative exercise or stress relief.
Application	<p>The Theraflex™ Putty can be used for rehab as well as for sports training or stress relief.</p> <ul style="list-style-type: none"> <li>• Rehabilitate post-operative / post-traumatic hands and arms</li> <li>• Restore function in hands, arms and shoulder following stroke</li> <li>• Reduce pain and stiffness in arthritic fingers, hands and wrists</li> <li>• Relax tense, sore muscles and reduce stress</li> <li>• Strengthen hands and arms</li> </ul>

### Technical specifications

Productgroup	Hand Therapy Products
Composition	<p>Silicone Elastomer          Non-toxic, greaseless, hypo-allergenic and latex-free          Oderless, stainless and will not leave any residue on patient's hands</p>

### Storage conditions

Duration of storage	Valid for 3 years from date of manufacture
Storage	Store in cool, dry place, avoid excessive heat (40° C.)


Product overview					
Art.Nr	Resistance	Color	Package	Weigt Package	Dimensions Ø x H
TP3PCT	Extra Soft	Tan	1 piece / cansiter	85 gr	Ø 8,2 x 4 cm
TP3PCY	Soft	Yellow	1 piece / cansiter	85 gr	Ø 8,2 x 4 cm
TP3PCR	Medium	Red	1 piece / cansiter	85 gr	Ø 8,2 x 4 cm
TP3PCG	Firm	Green	1 piece / cansiter	85 gr	Ø 8,2 x 4 cm
TP3PCB	Extra Firm	Blue	1 piece / cansiter	85 gr	Ø 8,2 x 4 cm
TP4PCT	Extra Soft	Tan	1 piece / cansiter	110 gr	Ø 11,7 x 4 cm
TP4PCY	Soft	Yellow	1 piece / cansiter	110 gr	Ø 11,7 x 4 cm
TP4PCR	Medium	Red	1 piece / cansiter	110 gr	Ø 11,7 x 4 cm
TP4PCG	Firm	Green	1 piece / cansiter	110 gr	Ø 11,7 x 4 cm
TP4PCB	Extra Firm	Blue	1 piece / cansiter	110 gr	Ø 11,7 x 4 cm
TP6PCT	Extra Soft	Tan	1 piece / cansiter	170 gr	Ø 11,7 x 4 cm
TP6PCY	Soft	Yellow	1 piece / cansiter	170 gr	Ø 11,7 x 4 cm
TP6PCR	Medium	Red	1 piece / cansiter	170 gr	Ø 11,7 x 4 cm
TP6PCG	Firm	Green	1 piece / cansiter	170 gr	Ø 11,7 x 4 cm
TP6PCB	Extra Firm	Blue	1 piece / cansiter	170 gr	Ø 11,7 x 4 cm
TPBLK1T	Extra Soft	Tan	1 piece / cansiter	450 gr	Ø 11,7 x 6,7 cm
TPBLK1Y	Soft	Yellow	1 piece / cansiter	450 gr	Ø 11,7 x 6,7 cm
TPBLK1R	Medium	Red	1 piece / cansiter	450 gr	Ø 11,7 x 6,7 cm
TPBLK1G	Firm	Green	1 piece / cansiter	450 gr	Ø 11,7 x 6,7 cm
TPBLK1B	Extra Firm	Blue	1 piece / cansiter	450 gr	Ø 11,7 x 6,7 cm


## Application

### Safety Precautions

Use the Theraflex™ Putty only after consulting a trained licensed healthcare professional.

## Exercises



Special items	
Customs code	9021.9090
Country of origin	U.S.A. MSD Europe is the authorised European Representative.
	This product has been manufactured according to CE-mark regulations.
Revision number	PI-Theraflex Putty (04.2008)
<b>MSD Europe bvba</b> Neringstraat 7 1840 Londerzeel Belgium Tel. +32 52 311 452 Fax. +32 52 311 450 info@msd-europe.com www.msd-europe.com	