

Touch-Test Discriminator

Measuring of cutaneous sensory



Characteristics

Function	<p>The Touch-Test Discriminator tests static and moving two-point discrimination of fingers, toes or sensory flaps. Consists of two sturdy, rotating, plastic disks that are joined together. The rounded tips are spaced at standard testing intervals.</p> <p>Ideal for testing following nerve repair, grafts and innervated tissue transfer for desensitization, or to determine the level of impairment.</p>
----------	---

Technical specifications

Productgroup	Evaluation line Instruments
Composition	plastic


Art.Nr	Description	Package	Weight Package	Dimensions L x W x H
12-1495	Touch-Test Discriminator	1 piece / sac	0,048 kg	

Storage conditions

Duration of storage	No limitation
Storage	No special requirements
Cleaning	Rince easily with water and soap

Application

Safety Precautions	Use the Touch-Test Discriminator only after consulting a trained licensed healthcare professional.
--------------------	--

Special items	
Customs code	9031.8098
Country of origin	U.S.A. MSD Europe is the authorised European Representative.
	This product has been manufactured according to CE-mark regulations.
Revision number	PI-Touch-Test Discriminator (04.2008)
MSD Europe bvba Neringstraat 7 1840 Londerzeel Belgium Tel. +32 52 311 452 Fax. +32 52 311 450 info@msd-europe.com www.msd-europe.com	