



MSD Instant Cold Pack

Instant cold therapy treatment

Characteristics

Function	<p>The MSD Instant Cold Pack is a disposable pack, one-time use only. This pack is used for immediate cold therapy.</p> <p>By finger pressure the pack will activate. The pack supplies instant cold to reduce heat and relieve pain and swelling, it stops little haemorrhages.</p> <p>The MSD Instant Cold Pack can be used in any first aid occasion, health care, daily use, travel and so on.</p>
----------	--

Technical specifications

Productgroup	Hot / Cold Therapy Products
Composition	Ammonium Nitrate and Water

Product overview


Art.Nr	Description	Package	Weight Package	Dimensions L x W x H
MSD-160	MSD Instant Cold Pack	10 pieces / box	2,7 kg	18 x 24 x 15 cm

Storage conditions

Duration of storage	Valid for 3 years from date of manufacture
Storage	Store in cool, dry place, avoid excessive heat (40° C.)
Cleaning	Rince easily with water and soap

Application

Safety Precautions	Use the MSD Instant Cold Pack only after consulting a trained licensed healthcare professional.
--------------------	---

Special items	
Customs code	3005.9099
Country of origin	China MSD Europe is the authorised European Representative.
	This product has been manufactured according to CE-mark regulations.
Revision number	PI-MSD Instant Cold Pack (04.2008)
MSD Europe bvba Neringstraat 7 1840 Londerzeel Belgium Tel. +32 52 311 452 Fax. +32 52 311 450 info@msd-europe.com www.msd-europe.com	